Falling Ball Viscometer

... Newtonian measurements made simple and easy!



What's Included?

Instrument Set of six (6) balls Temperature Probe Carrying Case

Optional Accessories

Temperature Bath (p33-35) Viscosity Standards (p52) Special Temperature Probes

Applications

Beverages

Coatings

Cosmetics

Detergents

Food

Paint

Petroleum Products

Pharmaceuticals

Polymers

Soap

How It Works

The Falling Ball Viscometer is based on the measuring principle by Höppler for simple but precise dynamic viscosity measurement of transparent Newtonian fluids. The basic concept is to measure the elapsed time required for the ball to fall under gravity through a sample-filled tube inclined at an angle*. The tube is mounted on a pivot bearing which quickly allows rotation of the tube 180 degrees, thereby allowing a repeat test to run immediately. Three measurements are taken and the average time it takes for the ball to fall is the result. A conversion formula turns the time reading into a final viscosity value.

The Falling Ball Viscometer is used for quality control in various industries as well as in academic institutions to illustrate scientific method. The ease of use and straightforward method for recording time measurements ensures meaningful test results.

*Model KF30 has a fixed angle of 80 degrees; Model KF40 can be angled at 50, 60, 70 and 80 degrees.

SPECIFICATIONS	
Viscosity Range:	0.5 mPa•s (cP) to 70,000 mPa•s (cP)
Accuracy:	0.5% - 2.0% depending on choice of ball
Ball set Material of Construction:	
	Balls 1 and 2: Boron Silicate Glass
	Balls 3 and 4: Nickel-iron
	Balls 5 and 6: Steel
Ball Diameter:	11.0 mm to 15.81 mm
Fall Time of Ball in Measurement: 30 to 300 seconds**	
Length of Measurement Zone in the Tube: 100 mm	
Operating Temperature Range: -5°C to +150°C	
Sample Tube Volume: 40mL	
Viscometer Dimensions: 180 x 220 x 330 mm	

^{**}Falling times greater than 300 seconds allow measurement of liquids above 70,000 mPa*s (cP)



Ball Set with Case

The Falling Ball Viscometer comes complete with a set of six (6) balls. See Specifications for material construction.



KF40 with Bath

Use with a Brookfield Circulating Bath permits rapid temperature control of sample for more accurate and repeatable results.