

培養箱 / 台灣製



These incubators are designed for drying and staining of slides tissue culture work, incubation of antibody tests, curing, pharmaceutical stability test, microbiological determinations, crystallization studies, incubation of hydroxysteroids, and similar application.

Specifications:

MODEL	MDI - 23	MDI - 80	MDI - 150
Temperature			
Working range	R.T. +5℃~85℃	R.T. +5℃~80℃	R.T. +5℃~80℃
Stability	±0.1℃	±0.1℃	±0.1℃
Uniformity	±0.5℃	±2℃	±2℃
Setting& readout	Digital	Digital	Digital
Controller	Microprocessor based	Microprocessor based	Microprocessor based
Heater	300 W	600 W	700 W
Shelf	1	2	2
Capacity	23 L	80 L	150 L
Dimension (cm)			
Chamber	W35 / D25 / H27	W40 / D40 / H50	W50 / D50 / H60
Overall	W45 / D44.5 / H47	W57 / D53 / H83	W67 / D63 / H92.5
Power source	110/220V 60/50Hz 5/3A	110/220V 60/50Hz 6/3A	110/220V 60/50Hz 7/4A
Weight	22 kg	70 kg	90 kg

明頂科技有限公司

<https://www.mding.com.tw>

TEL.07-552-6160 FAX.07-586-6000